



PRODUCT CATALOGUE

YOUR
PARTNER
IN DIAGNOSTIC
EXCELLENCE

Technology and Innovation to Improve the Quality of Life

Qingdao Hightop Biotech Co., Ltd.

CONTENTS

01

HIGHTOP Milestones

01

02

Certifications

02

03

Hematology

03

04

Urinalysis

04

05

Biochemistry

05-06

06

Chemiluminescence
Immunoassay (CLIA)

07-08

07

ELISA

09-10

08

Point of Care (POCT)

11-13

09

Rapid Test

14-24

9.1

Fertility Test

9.2

Infectious Disease Tests

9.3

Drug of Abuse Tests

9.4

Cardiac Marker Tests

9.5

Tumor Marker Tests

9.6

Alcohol Tests

9.7


Inflammation Tests

9.8

Other Rapid Tests

9.9

DOT Kits



ABOUT HIGHTOP

2009

HIGHTOP established

Founded in Qingdao 2009, HIGHTOP has been focused on innovating, manufacturing, and marketing of medical devices and in vitro diagnostic reagents ever since.



1996

Kanghua established

HIGHTOP is the subsidiary of Shandong Kanghua Biotech Co., Ltd. Which founded in 1996, and is welcoming its 30th anniversary. Together with Kanghua, we are now owning over 600 R&D engineers and have developed a wide range of diagnostic product line that are using in over 200 countries and territories, including Hematology, Biochemistry, Chemiluminescence, Immunofluorescence, Urinalysis, Colloidal Gold and DOT rapid test, ELISA, etc.

600+

R&D Engineer

200+

Countries and Regions

We value our customers as life, devote ourselves to supply customers the most competitive products and considerate professional service. At HIGHTOP, Where Quality Meets User Satisfaction, we are dedicated to become Your Partner in Diagnostic Excellence!

HIGHTOP MILESTONES

- 1996** Kanghua established
- 2002** Built the first independent R&D center
- 2004** Our company passed ISO9001 and ISO13485 certifications
- 2004** Self-developed hCG & LH Test Strips were certified by NMPA
- 2009** First NMPA certified Colloidal Gold MP-IgM Antibody Test Kit
- 2010** Built the Fluorescence Immunoassay R&D team
- 2010** TORCH panels in ELISA and Colloidal Gold were certified by NMPA
- 2013** Established CLIA R&D team
- 2014** Built one of the biggest IVD industrial park in China
- 2016** Successfully obtained TÜV SÜD ISO13485 certification and EU CE certification of 39 products, officially expanding into international markets
- 2017** First VEGF reagent on CLIA
- 2019** Self-developed & self-produced for core raw materials
- 2019** HIGHTOP hCG Rapid Test has become the company's first FDA approved product
- 2023** HIGHTOP acquired MDSAP certification
- 2024** First NMPA certified colloidal gold respiratory 7 in 1 rapid test panel
- 2025** Ultra-high Throughput CLIA Analyzer Aurora-6000i launched
- 2025** HIGHTOP hCG & LH Rapid Tests has become the company's first batch CE-IVDR certificated products

CERTIFICATIONS



800+
Patents

370+
CE

4
FDA

80+
ANVISA



Hematology CE

3-part differentiation

Support 21 items, 3 histograms;
Stores 50,000 results

Highly Efficient

- 2 Reagents, sampling volume $\leq 20\mu\text{L}$
- 60 samples/hour
- Average 59 seconds/sample

Compact Yet Smart

- Support connection to LIS
- Built-in thermal printer and external printer both available
- Size: 390×285×390 mm; Weight $\leq 18\text{kg}$
- Support multilanguage system

HT-300



HT-300VET



Support 13 species of animals: Cat, Dog, Pig, Rabbit, etc.

5-part differentiation

Support 25 items, 2 histograms, and 2 scatterplots;
Stores 50,000 results

Compact Yet Smart

- Support connection to LIS
 - Support External printer
 - Size: 430×330×510 mm
- Weight $\leq 25\text{kg}$

Highly Efficient

- 3 Reagents, sampling volume $\leq 28\mu\text{L}$
- 55 Samples/hour

HT-5000



Veterinary use available:
Support 13 species of animals

Lyse (HT-HE) : 100mL/bottle, 200mL/bottle, **500mL/bottle**

LySe (B-DIFF): 500mL/bottle, **1L/bottle**, 2L/bottle


Diluent: 5L/bottle, 10L/bottle, **20L/bottle**

EZ Cleaner (optional): 25mL/bottle

Probe Cleaning (optional): 25mL/bottle



Urinalysis CE NMPA



The Urinalysis Reagent Strip can test for 14 parameters:

Uncut sheet available of any combination of following parameters

Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocytes, Ascorbic Acid, in urine specimen, Calcium, Creatinine and Microalbuminuria.

compatible QC reagent available on-demand



HT-115A

Throughput	514 tests/hour
Strips	4-14 parameters
Data Storage	5,000



HT-111A

Detection Speed	60/120 reports/hour
Strips	4-14 parameters
Data Storage	5,000



HT-1420A

Detection Speed	240 tests/hour
Strips	10-14 parameters
Data Storage	10,000

Standard Size: 20T/50T/100T

Cat. No.	Product Description
H411	URS-1G(GLUCOSE)
H412	URS-1H(pH)
H413	URS-1K(KETONE)
H414	URS-1A(ASCORBIC)
H415	URS-1L(LEUKOCYTES)
H416	URS-1M(MICROALBUMINURIA)
H417	URS-1Bi(BILIRUBIN)
H418	URS-1Bl(BLOOD)
H419	URS-1S(SPECIFICGRAVITY)
H4110	URS-1P(PROTEIN)
H4111	URS-1N(NITRITE)
H4112	URS-1U(UROBILINOGEN)
H4113	URS-1Cr(CREATININE)
H4114	URS-1Ca(CALCIUM)
H421	URS-2K(GLUCOSE/KETONE)
H422	URS-2P(GLUCOSE/PROTEIN)

Cat. No.	Product Description
H423	URS-2N(LEU/NIT)
H424	URS-2M(CR/MALB)
H425	URS-2C(CR/Ca)
H426	URS-2B(URO/BIL)
H431	URS-3K(GLUCOSE/PROTEIN/KET)
H432	URS-3H(GLUCOSE/PROTEIN/pH)
H441	URS-4 (GLUCOSE/PROTEIN/pH/SG)
H442	URS-4B (GLUCOSE/PROTEIN/pH/BLOOD)
H451	URS-5(GLU/PRO/pH/BLD/KET)
H452	URS-5H(GLU/PRO/KET/BLD/PH)
H491	URS-9(LEU/NIT/URO/PRO/pH/BLD/KET/BIL/GLU)
H4101	URS-10S(LEU/NIT/URO/PRO/pH/BLD/SG/KET/BIL/GLU)
H4102	URS-10C(LEU/NIT/URO/PRO/pH/BLD/KET/BIL/GLU/VC)
H4111	URS-11(URO/BLD/BIL/KET/LEU/GLU/PRO/pH/NIT/SG/VC)
H4141	URS-14(VC/LEU/GLU/MALB/BIL/pH/BLD/SG/KET/PRO/NIT /URO/CR/CA)

Biochemistry Analyzer CE

HTSH-2000

200 Tests/hour

Size 622×428×527mm

Weight 47.5KG



Reagent sample tray

40 Reagent Positions, 48 Sample Positions
Mixed display of all kinds of samples.
37°C ± 0.1°C, real-time monitor.

Reagent/sample filling mechanism

Special coating to avoid contamination.

Reaction tray

50 Reaction Positions
Automatically monitor on the remaining Circulating cooling refrigeration.

Stirring mechanism

Liquid level sensor function. Anti-collision function.
Reagent volume real-time detection.
Automatic cleaning function.

Automatic cleaning mechanism

Automatic eight sections washing of detecting cup.
Washing probe has the function of plugging and automatic cleaning.


HTSH-4000

• Test Speed	400 Tests/Hour
• Sample Position	96
• Reagent Position	80
• Reaction Position	90
• Reaction Volume	150~500 µL
• Dimension	930 x 700 x 1130mm (LxWxH)
• Weight	112 Kg



HTSH-8000

• Test Speed	800 Tests/Hour
• Sample Position	144 (93+15+18+15+3)
• Reagent Position	100
• Reaction Position	165
• Reaction Volume	90~450 µL
• Dimension	1568 x 892 x 1120mm (LxWxH)
• Weight	400 Kg



Biochemistry Reagents NMPA

	Product Description	Qualification
Liver Function Tests	ALB Albumin	CE
	TP Total Protein	CE
	ALP Alkaline Phosphatase	CE
	ALT Alanine Aminotransferase	CE
	AST Aspartate Aminotransferase	CE
	γ-GT Gamma-Glutamyl Transpeptidase	CE
	DB Direct Bilirubin	CE
	TB Total Bilirubin	CE
	CHE Choline Esterase	
	AFU α-L-Fucosidase	
	TBA Total Bile Acid	
	ADA Adenosine Deaminase	
	5'-NT 5'-Nucleotidase	
	mAST Mitochondrial Aspartate Aminotransferase	
	CG Cholyglycine	
	PA Prealbumin	
	MAO Monoamine Oxidase	
	GR Glutathione Reductase	
Blood-lipids Tests	CHOL Total Cholesterol	CE
	TG Triglycerides	CE
	HDL-C High Density Lipoprotein Cholesterol	CE
	LDL-C Low Density Lipoprotein Cholesterol	CE
	NEFA Nonesterified Fatty Acid	
	Lp(a) Lipoprotein(a)	
	ApoA1 Apolipoprotein A1	
	ApoB Apolipoprotein B	
	ApoE Apolipoprotein E	
	HCY Homocysteine	
	sdLDL-C Small Dense Low-Density Lipoprotein Cholesterol	
	Lp-PLA2 Lipoprotein-associated Phospholipase A2	
Glucose Metabolism Tests	GLU Glucose	CE
	FMN Serum Fructosamine	
	HbA1c Glycosylated Hemoglobin	
	GA Glycated Albumin	
	ADPN Adiponectin	
	β-HB β-hydroxybutyrate	
Renal Function Tests	Urea Urea	CE
	UA Uric Acid	CE
	Crea Creatinine	CE
	CYS-C Cystatin C	
	β ₂ -MG β ₂ -Microglobulin	
	α ₁ -MG α ₁ -microglobulin	
	NGAL Neutrophil Gelatinase-Associated Lipocalin	
	NAG N-acetyl -β-D-glucosaminidase	
	MALB Urine Microalbumin	
	RBP Retinol-binding Protein	
Pancreatic Gland Tests	α-AMY α-Amylase	
	LPS Lipase	
Myocardial Function Tests	CK Creatine Kinase	
	CK-MB Creatine Kinase Isoenzyme	
	LDH Lactate Dehydrogenase	
	α-HBDH α-hydroxybutyrate Dehydrogenase	
	LDH1 Lactic Acid Dehydrogenase Isozyme 1	
	Myo Myoglobin	
	cTnI Cardiac Troponin I	
	H-FABP Heart-type Fatty Acid-Binding Protein	
	IMA Ischemia Modified Albumin	
	ACE Angiotensin Converting Enzyme	
Electrolyte & Blood Gase Tests	Ca Calcium	
	P Phosphorus	
	CO ₂ Carbon Dioxide	
	Cl Chloride	
	Mg Magnesium	
	Fe Iron	
	Cu Cuprum	
	Zn Zinc	
	IgA Immunoglobulin A	
	IgG Immunoglobulin G	
Specialty Protein Tests	IgM Immunoglobulin M	
	C3 Complement C3	
	C4 Complement C4	
	Fg Fibrinogen	
	AAG Alpha1-acidic Glycoprotein	
	HPT Haptoglobin	
	TRF Transferrin	
	AAT Alpha1-antitrypsin	
	IgG4 Immunoglobulin G4	
	C1q Complement C1q	
	ASO Antistreptolysin"O"	
	RF Rheumatic Factor	
	CRP C-reaction Protein	
	hs-CRP High-Sensitivity C-reactive Protein	
	Anti-CCP Anti-Cyclic Citrullinated Peptide	
	PCT Procalcitonin	
	SAA Serum Amyloid A	
	MMP3 Metallo Matrix Protease 3	
	KL-6 Krebs von den Lungen-6 antigen	
Others Tests	PGI Pepsinogen I	
	PGII Pepsinogen II	
	SOD Superoxide Dismutase	
	SA Sialic Acid	
	CMCAL Multiple biochemical calibrators(Normal value)	CE
Calibration and Control	CMCON1 Biochemical multiple quality control material(Level 1)	CE
	CMCON2 Biochemical multiple quality control material(Level 2)	CE

Chemiluminescence Immunoassay (CLIA)

Aurora-1000i

225 Tests/hour

Size 1060×590×720mm

Weight 105kg



Acridinium Ester labeled direct chemiluminescence

700	700 reaction vessel capacity for continuous loading	20	20 refrigerated reagent positions (2-8 °C)
60	60 sample positions with STAT function	86	More than 86 reagents available

Aurora-500i

• Detection speed	72 tests/hour
• Reagent position	32
• Sample position	10
• Simultaneously analyzed reagents	8
• Reaction cup capacity	400
• Time to first result	11 minutes



Aurora-2000i

• Detection speed	200 tests/hour
• Reagent position	75
• Sample position	125
• Incubation position	136
• Simultaneously analyzed reagents	25
• Reaction cup capacity	1000



Aurora-6000i

• Detection speed	600 tests/hour, up to 2400 tests/hour
• Reagent position	50
• Sample position	310
• Incubation position	288
• Simultaneously analyzed reagents	50
• Reaction cup capacity	3000



Chemiluminescence Immunoassay (CLIA)

Kit Size 50T/100T

	Product Description	Qualification
Prenatal Screening	Pregnancy AFP	CE NMPA
	PAPP-A	CE NMPA
	Free β -HCG	CE NMPA
	uE3	CE NMPA
Thyroid	T3	CE NMPA
	T4	CE NMPA
	FT3	CE NMPA
	FT4	CE NMPA
	TSH	CE NMPA
	TG-Ab	CE NMPA
	TPO-Ab	CE NMPA
	TG	CE NMPA
Tumor Marker	Ferritin	CE NMPA
	VEGF	CE NMPA
	S100	CE NMPA
	PG-I	CE NMPA
	PG-II	CE NMPA
	AFP	CE NMPA
	CEA	CE NMPA
	G-17	CE NMPA
	CA15-3	CE NMPA
	CA125	CE NMPA
	CA19-9	CE NMPA
	PSA	NMPA
	f-PSA	NMPA
	CA50	CE NMPA
	CA72-4	CE NMPA
	NSE	CE NMPA
	β 2-MG	CE NMPA
	CA242	CE NMPA
	SCC	CE NMPA
	proGRP	CE NMPA
	CYFRA21-1	CE NMPA
	HE4	CE NMPA
	TPS	NMPA
	S100 β	NMPA
	PACP	NMPA
Fertility	LH	CE NMPA
	FSH	CE NMPA
	PRL	CE NMPA
	E2	CE NMPA
	Progesterone	CE NMPA
	Testosterone	CE NMPA
	HCG	CE NMPA
	Total β -HCG	CE NMPA
Diabetes	INS	NMPA
	C-P	NMPA

	Product Description	Qualification
Bone Metabolism	25-OH Vitamin D	CE NMPA
	Intact PTH	CE NMPA
Anemia	Folate	CE NMPA
	Vitamin B12	CE NMPA
Inflammation	sTfR	NMPA
	PCT	CE NMPA
	IL-6	CE NMPA
	SAA	CE NMPA
	CRP	CE NMPA
	Lp-PLA2	NMPA
	hs-cTnI	NMPA
	NT-proBNP	CE NMPA
Cardiac	MYO	CE NMPA
	h-FABP	CE NMPA
	cTnI	CE NMPA
	CK-MB	CE NMPA
	MPO	NMPA
	D-Dimer	CE NMPA
	BNP	NMPA
	sST2	NMPA
	MP IgG	NMPA
	MP IgM	NMPA
Respiratory Diseases	Cpn IgG	NMPA
	Cpn IgM	NMPA
	RSV IgM	
	Adv IgM	
	FLU-A IgM	
	FLU-B IgM	
	CBV IgM	
	HIV Ag/Ab	NMPA
Infectious Disease	HCV Ab	NMPA
	TP Ab	NMPA
	HBsAg	NMPA
	HBsAb	NMPA
	HBeAg	NMPA
	HBeAb	NMPA
	HBcAb	NMPA
	CG	NMPA
Liver Fibrosis	HA	NMPA
	LN	NMPA
	CIV	NMPA
	PIIIN-P	NMPA
	FN	NMPA
Renal	PLA2R IgG	NMPA

	Cat. No.	Product Description	Qualification
Blood Screening Tests	H3113	HIV 1/2 Antibody(3rd) Test	
	H3114	HIV 1/2 Antibody/Antigen(4th) Test	
	H312	HBsAg Hepatitis B Surface Antigen Test	
	H313	HCV Hepatitis C Virus Antibody Test	
	H314	Syphilis Antibody Test	CE
	H315	HBsAb Hepatitis B Surface Antibody Test	
	H316	HBeAg Hepatitis B Envelope Antigen Test	
	H317	HBeAb Hepatitis B Envelope Antibody Test	
	H318	HBcAb Hepatitis B Core Antibody Test	
Quantitative Tests	H321	AFP Alpha Fetoprotein Antigen Test	CE NMPA
	H322	CEA Carcinoma Embryonic Antigen Test	CE NMPA
	H323	T3 3,5,3'-Triiodo-Thyronine Test	CE
	H324	T4 Thyroxine Test	CE
	H325	FT3 Free 3,5,3'-Triiodo-Thyronine Test	
	H326	FT4 Free Thyroxine Test	
	H327	TSH Thyroid Stimulating Hormone Test	CE
	H328	LH Luteinizing Hormone Test	CE
	H329	FSH Follicle Stimulating Hormone Test	CE
	H3210	SF Serum Ferritin Test	CE
	H3211	PSA Prostate Specific Antigen Test	
	H3212	VEGF Vascular Endothelial Growth Factor Test	CE
	H3213	Total IgE Antibody Test	CE
	H3214	β 2 Glycoprotein Antibody Test	CE
	H3215	AMH Anti-Mullerian Hormone Test	CE NMPA
	H3216	PRL Prolactin Test	
Prenatal and Postnatal Care Tests	H331M	TOXO Toxoplasma Ab IgM Test	NMPA
	H331G	TOXO Toxoplasma Ab IgG Test	NMPA
	H332M	RV Rubella Virus Ab IgM Test	NMPA
	H332G	RV Rubella Virus Ab IgG Test	NMPA
	H333M	CMV Cytomegalovirus Ab IgM Test	NMPA
	H333G	CMV Cytomegalovirus Ab IgG Test	NMPA
	H334M	HSV-I Herpes Simplex Virus I Ab IgM Test	CE NMPA
	H334G	HSV-I Herpes Simplex Virus I Ab IgG Test	CE NMPA
	H335M	HSV-II Herpes Simplex Virus II Ab IgM Test	CE NMPA
	H335G	HSV-II Herpes Simplex Virus II Ab IgG Test	CE NMPA
Hepatitis Virus Tests	H336M	HSV1&2 Herpes Simplex Virus (HSV1&2) Ab IgM Test	
	H336G	HSV1&2 Herpes Simplex Virus (HSV1&2) Ab IgG Test	
	H351M	HAV Hepatitis A Virus Ab IgM Test	CE NMPA
Other Tests	H352M	HEV Hepatitis E Virus Ab IgM Test	CE NMPA
	H361M	TB Tuberculosis Ab IgM Test	CE
	H361G	TB Tuberculosis Ab IgG Test	CE NMPA
	H362G	ds-DNA Anti-Double Stranded DNA Ab IgG Test	CE NMPA
	H363G	ANA Anti-Nuclear Ab IgG Test	CE NMPA
	H364	HP H.Pylori Antigen Test	CE
	H365G	Ct Chlamydia Trachomatis Ab IgG Test	
	H365M	Ct Chlamydia Trachomatis Ab IgM Test	
	H366	SARS-CoV-2 Neutralizing Antibody Test	CE

Immunological Infertility Tests

Cat. No.	Product Description	Qualification
H341A	ASA Anti-Spermatozoa Ab IgA Test	CE NMPA
H341M	ASA Anti-Spermatozoa Ab IgM Test	CE NMPA
H341G	ASA Anti-Spermatozoa Ab IgG Test	CE NMPA
H342A	EMAb Anti-Endometrial Ab IgA Test	CE
H342M	EMAb Anti-Endometrial Ab IgM Test	CE NMPA
H342G	EMAb Anti-Endometrial Ab IgG Test	CE NMPA
H343A	ACA Anticardiolipin Ab IgA Test	CE NMPA
H343M	ACA Anticardiolipin Ab IgM Test	CE NMPA
H343G	ACA Anticardiolipin Ab IgG Test	CE NMPA
H344A	AOA Anti-Ovarian Ab IgA Test	CE
H344M	AOA Anti-Ovarian Ab IgM Test	CE NMPA
H344G	AOA Anti-Ovarian Ab IgG Test	CE NMPA
H345A	HCG Anti-Human Chorionic Gonadotrophin Ab IgA Test	CE
H345M	HCG Anti-Human Chorionic Gonadotrophin Ab IgM Test	CE
H345G	HCG Anti-Human Chorionic Gonadotrophin Ab IgG Test	CE
H346A	ZP Anti-zona Pellucida Ab IgA Test	CE
H346M	ZP Anti-zona Pellucida Ab IgM Test	CE
H346G	ZP Anti-zona Pellucida Ab IgG Test	CE
H347A	TA Anti-Trophoblast Cell Membrane Ab IgA Test	CE
H347M	TA Anti-Trophoblast Cell Membrane Ab IgM Test	CE
H347G	TA Anti-Trophoblast Cell Membrane Ab IgG Test	CE
H348M	EHF Epidemic Hemorrhagic Fever Virus Ab IgM Test	CE NMPA



Point of Care (POCT)

Fluorescence Immunoassay System

RapidFF-1

Data Transmission	Supports multiple transmission mode
Size	112×65×40mm
Weight	0.2kg with battery
Test Speed	720 Tests/hour
Printer	Supports connection to mobile mini program or PC software for printing



HTY-100

Data Transmission	Connected to LIS by COM/USB
Screen	LCT Touch Screen
Printer	Internal Thermal Printer
Size	264×233×142mm
Weight	About 2.2kg
Test Speed	Fastest 3 second/test, routinely 3 min/test



HTY-1600

Data Transmission	Connected to LIS by COM/USB
Screen	19.68 in LCT Touch Screen 15inch
Printer	Internal Thermal Printer & External Printer
Size	430×414×506mm
Weight	About 27kg
Test Speed	200 Tests/hour



RapidPro6000

Test Speed	170 tests/hour
Sample Position	30
Reagent Position	5×50
Time to First Result	3 min
Size	709×747×758mm
Weight	85kg



Point of Care (POCT)

	Cat. No.	Diseases	Test Item	Specimen	Qualification	Kit Size
Inflammation	H611	Inflammation	CRP C-reactive Protein	WB/S/P/fingertip Blood	CE NMPA	25T
	H612		PCT Procalcitonin	WB/S/P	CE NMPA	25T
	H613		SAA Serum Amyloid A	WB/S/P/fingertip Blood	CE NMPA	25T
	H614		CRP/SAA Combo	WB/S/P/fingertip Blood	CE NMPA	25T
	H615		HBP Heparin-binding Protein	Whole Blood/Plasma	CE	25T
	H616		CCP Cyclic Citrullinated Peptide IgG	WB/S/P	CE NMPA	25T
	H617		RF/ASO/CRP Combo	Serum	CE NMPA	25T
	H618		RF Rheumatoid factor	WB/S/P	CE	25T
	H619		ASO Anti-streptolysin O	WB/S/P	CE	25T
	H6110		IL-6 Interleukin 6	WB/S/P	CE NMPA	25T
Cardiac Markers	H621	Myocardial Infarction	cTnI Cardiac Troponin I	WB/S/P	CE NMPA	25T
	H622		CK-MB Creatine Kinase Muscle B	WB/S/P	CE NMPA	25T
	H623		Myo Myoglobin	WB/S/P	CE NMPA	25T
	H626		H-FABP Heart-type Fatty Acid-Binding Protein	Serum/Plasma	CE NMPA	25T
	H627		cTnI/CK-MB/Myo Combo	WB/S/P	CE NMPA	25T
	H6212		cTnT Cardiac Troponin T	WB/S/P	CE NMPA	25T
	H628		cTnI/NT-proBNP/D-D Combo	Whole Blood/Plasma	CE NMPA	25T
	H6214		hs-cTnI High Sensitive Cardiac Troponin I	WB/S/P	NMPA	25T
	H6215		CPP Copeptin	WB/S/P		25T
	H6216		GDF-15 Growth differentiation factor-15	WB/S/P		25T
	H6211		Lp-PLA2 Lipoprotein-associated phospholipaseA2	WB/S/P	CE NMPA	25T
	H629		ST2 Growth STimulation expressed gene 2	WB/S/P	CE NMPA	25T
	H6210		NT-proBNP/ST2 Combo	WB/S/P	CE NMPA	25T
	H624		D-D D-dimer	Whole Blood/Plasma	CE NMPA	25T
	H625		NT-proBNP N-terminal brain natriuretic peptide	WB/S/P	CE NMPA	25T
	H6213		BNP B-type brain natriuretic peptide	Whole Blood/Plasma	CE NMPA	25T
Tumor Markers	H631	Liver Cancer	AFP Alpha-fetoprotein	Serum/Plasma	CE NMPA	25T
	H632	Prostate Cancer	t-PSA Prostate Specific antigen	WB/S/P	CE NMPA	25T
	H633	Cancer	CEA Carcinoembryonic antigen	Serum/Plasma	CE NMPA	25T
	H634	Gastric Cancer	PG I Pepsinogen I	Serum	CE NMPA	25T
	H635		PG II Pepsinogen II	Serum	CE NMPA	25T
	H6310		PGI+PGII Pepsinogen I/Pepsinogen II	Serum		25T
	H636	Colorectal Cancer	FOB Fecal Occult Blood	Feces	CE	25T
	H637	Cerebral injury	S100 β Central nervous system specific protein	WB/S/P	CE NMPA	25T
	H638	Vascular hyperplasia	VEGF Vascular Endothelial Growth Factor	WB/S/P	CE NMPA	25T
	H639	Atrophic gastritis	G-17 Gastrin 17	WB/S/P	CE NMPA	25T
	H6311	Ovarian Cancer	CA125 Carbohydrate Antigen 125	WB/S/P		25T
	H6312	Breast Cancer	CA153 Carbohydrate Antigen 15-3	WB/S/P		25T
	H6313	Pancreatic Cancer	CA199 Carbohydrate Antigen 199	WB/S/P		25T

	Cat. No.	Diseases	Test Item	Specimen	Qualification	Kit Size	
Infectious Diseases	H651	Respiratory Disease	MP IgM Mycoplasma pneumoniae IgM	WB/S/P	CE NMPA	25T	
	H652		MP IgG Mycoplasma pneumoniae IgG	WB/S/P	CE NMPA	25T	
	H653		Cpn IgM Chlamydia pneumoniae IgM	WB/S/P	NMPA	25T	
	H654		SARS-CoV-2 Neutralizing Antibody	WB/S/P	CE	25T	
	H655		Cpn IgG Chlamydia pneumoniae IgG	WB/S/P	NMPA	25T	
	H656		MP Mycoplasma Pneumoniae	Oropharyngeal swab	CE	25T	
	H657		Cpn Chlamydia pneumoniae	Oropharyngeal swab		25T	
	H658		Adv Adenovirus Antigen	Nasopharynx swab	CE	25T	
	H659		RSV Respiratory Syncytial Virus Antigen	Nasopharynx swab	CE	25T	
	H6510		Strep A Antigen	Oropharyngeal swab	CE	25T	
	H6511		Influenza A+B Antigen	Nasopharynx swab	CE	25T	
	H6512	Gastric Disease	H.Pylori Antibody	WB/S/P	CE	25T	
	H6513	Acute Febrile Disease	Dengue NS1 Antigen	WB/S/P	CE	25T	
	H6514		Dengue IgM/IgG Antibody	WB/S/P	CE	25T	
Hormone Markers	H661	Pregnancy	Total β-HCG	Serum	CE NMPA	25T	
	H667		FSH follicle-stimulating hormone	WB/S/P	CE	25T	
	H668		LH luteinizing hormone	WB/S/P	CE NMPA	25T	
	H669		HCG Human Chorionic Gonadotropin	Serum	CE	25T	
	H6610		AMH Anti-Müllerian Hormone	WB/S/P		25T	
	H6611		Prog Progesterone	WB/S/P		25T	
	H6612	Male Gonadal Abnormalities	Testo Testosterone	WB/S/P		25T	
	H6613	Growth Hormone	IGF-1 Insulin-Like Growth Factor-1	WB/S/P		25T	
	H6614		IGFBP-3 Insulin-like Growth Factor Binding Protein-3	WB/S/P		25T	
	H6615		IGF-1/IGFBP-3 Insulin-like Growth Factor-1	WB/S/P		25T	
			/Insulin-like Growth Factor Binding Protein-3				
	H662	Thyroid Function	T3 Triiodothyronine	Serum	CE	25T	
	H663		T4 Thyroxine	Serum	CE	25T	
	H664		TSH Thyrotrophic Hormone	Serum	CE	25T	
	H665		FT4 Free thyroxine	Serum	CE	25T	
	H666		FT3 Free triiodothyronine	Serum	CE	25T	
Nutritional Metabolism Markers	H671	25-OH-D Deficiency	25-OH-D 25-Hydroxy Vitamin D	WB/S/P	CE NMPA	25T	
	H681	Iron-deficiency Anemia	Ferr Ferritin	WB/S/P	CE NMPA	25T	
	H682	Folate Deficiency	FA Folic Acid	WB/S/P	CE NMPA	25T	
	H683	Vitamin B12 Deficiency	VB12 Vitamin B12	WB/S/P	CE NMPA	25T	
Diabetic & Renal Injury Markers	H641	Renal Injury	Cys C Cystatin C	Serum/Plasma	CE NMPA	25T	
	H642		β2-MG Beta 2-Microglobulin	Serum/Plasma	CE NMPA	25T	
	H643		NGAL Neutrophil Gelatinase Associated	Serum	CE NMPA	25T	
	H644		MAU Microalbumin	Urine	CE NMPA	25T	
	H645	Diabetes Mellitus	HbA1C Hemoglobin A1c	Whole Blood	CE NMPA	25T	
	H646		ADPN Adiponectin	WB/S/P		25T	
Others	H691	Digestive Disease	HB+TF Hemoglobin+Transferrin Combo	Feces	CE NMPA	25T	

Rapid Test

Fertility Test

hCG Pregnancy Test

The hCG Rapid Test is a rapid chromatographic immunoassay, used for early detecting human chorionic gonadotropin (hCG) in human urine and/or serum to indicate pregnancy.

Features

- Serum/urine hCG test in one strip/cassette
- Cut-off available: 10mIU/mL, 20mIU/mL, 25mIU/mL
- CE certified by European Notified Body TUV (CE0123)



Strip



Cassette



Midstream

	Cat. No.	Product Description	Qualification	Specimen	Cut-Off	Format	Kit Size
hCG Pregnancy Test	H111-11U	hCG Pregnancy Rapid Test	CE ₀₁₂₃ FDA NMPA	Urine	10mIU/mL, 25mIU/mL	Strip	50T
	H111-10U		CE ₀₁₂₃ FDA NMPA	Urine	10mIU/mL, 25mIU/mL	Strip	1T
	H111-21U		CE ₀₁₂₃ FDA NMPA	Urine	10mIU/mL, 25mIU/mL	Cassette	25T
	H111-20U		CE ₀₁₂₃ FDA NMPA	Urine	10mIU/mL, 25mIU/mL	Cassette	1T
	H111-31U		CE ₀₁₂₃ FDA NMPA	Urine	10mIU/mL, 25mIU/mL	Midstream	20T
	H111-30U		CE ₀₁₂₃ FDA NMPA	Urine	10mIU/mL, 25mIU/mL	Midstream	1T
	H111-11P		CE NMPA	Serum	10mIU/mL, 25mIU/mL	Strip	50T
	H111-21P		CE NMPA	Serum	10mIU/mL, 25mIU/mL	Cassette	25T
	H111-11B		CE NMPA	Urine/Serum	10mIU/mL, 25mIU/mL	Strip	50T
	H111-21B		CE NMPA	Urine/Serum	10mIU/mL, 25mIU/mL	Cassette	25T
	H111-32U	hCG Pregnancy Weekly Rapid Test	CE NMPA	Urine	-	Midstream	1T
LH Ovulation Test	H112-11U	LH Ovulation Rapid Test	CE ₀₁₂₃ FDA NMPA	Urine	25mIU/mL	Strip	50T
	H112-21U		CE ₀₁₂₃ FDA NMPA	Urine	25mIU/mL	Cassette	25T
	H112-31U		CE ₀₁₂₃ FDA NMPA	Urine	25mIU/mL	Midstream	20T
	H112-30U		CE ₀₁₂₃ FDA NMPA	Urine	25mIU/mL	Midstream	1T
FSH Menopause Test	H113-11U	FSH Menopause Rapid Test	CE	Urine	25mIU/mL	Strip	50T
	H113-21U		CE	Urine	25mIU/mL	Cassette	25T

Rapid Test

Infectious Disease Tests

	Cat. No.	Product Description	Qualification	Specimen	Format	Kit Size
Hepatitis Virus Tests	H121-21WG	HAV Hepatitis A Virus Ab IgG Test		WB/S/P	Cassette	25T
	H121-11WM	HAV Hepatitis A Virus IgM Ab Rapid Test	CE NMPA	WB/S/P	Strip	50T
	H121-21WM		CE NMPA	WB/S/P	Cassette	25T
	H121-21W	HAV Hepatitis A Virus IgG/IgM Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H121-41W		CE	WB/S/P	Panel	25T
	H122-41W	HBV Hepatitis B Virus 5-in-1 Rapid Test	NMPA	WB/S/P	Panel	20T
	H123-11W	HBsAg Hepatitis B Surface Ag Rapid Test	NMPA	WB/S/P	Strip	50T
	H123-21W		NMPA	WB/S/P	Cassette	25T
	H124-11W	HBsAb Hepatitis B Surface Ab Rapid Test		WB/S/P	Strip	50T
	H124-21W			WB/S/P	Cassette	25T
	H125-11W	HBeAg Hepatitis B Envelope Ag Rapid Test		WB/S/P	Strip	50T
	H125-21W			WB/S/P	Cassette	25T
	H126-11W	HBeAb Hepatitis B Envelope Ab Rapid Test		WB/S/P	Strip	50T
	H126-21W			WB/S/P	Cassette	25T
	H127-11W	HBcAb Hepatitis B Core Ab Rapid Test		WB/S/P	Strip	50T
	H127-21W			WB/S/P	Cassette	25T
	H128-11W	HCV Hepatitis C Virus Ab Rapid Test	NMPA	WB/S/P	Strip	50T
	H128-21W		NMPA	WB/S/P	Cassette	25T
	H129-11W	HEV Hepatitis E Virus Ab IgM Rapid Test	CE NMPA	WB/S/P	Strip	50T
	H129-21W		CE NMPA	WB/S/P	Cassette	25T
	H1210-41W	HBsAg/HCV Rapid Test		WB/S/P	Panel	20T
Sexually Transmitted Diseases Tests	H131-11W	HIV 1/2 Ab Rapid Test	NMPA	WB/S/P	Strip	50T
	H131-21W		NMPA	WB/S/P	Cassette	25T
	H131-11U			Urine	Strip	50T
	H131-21U			Urine	Cassette	25T
	H131-11O			Oral swab	Strip	50T
	H131-21O			Oral swab	Cassette	25T
	H131-11WT	HIV 1/2 Ab Rapid Test (Tri-line)		WB/S/P	Strip	50T
	H131-21WT			WB/S/P	Cassette	25T
	H132-11W	HIV Ag/Ab Rapid Test		WB/S/P	Strip	50T
	H132-21W			WB/S/P	Cassette	25T
	H1312-21W	HIV 1.2.0 Ab Rapid Test (Tri-line)		WB/S/P	Cassette	25T
	H133-11W	Syphilis (TP) Ab Rapid Test	CE NMPA	WB/S/P	Strip	50T
	H133-21W		CE NMPA	WB/S/P	Cassette	25T
	H134-21W	HIV/Syphilis Rapid Test (Tri-line)		WB/S/P	Cassette	25T
	H135-41W	HIV/HBsAg/HCV 3-in-1 Rapid Test		WB/S/P	Panel	20T
	H137-41W	HIV/HBsAg/HCV/Syphilis 4-in-1 Rapid Test		WB/S/P	Panel	20T
	H138-11S	Gonorrhoeae Ag Rapid Test	CE NMPA	Swab	Strip	50T
	H138-21S		CE NMPA	Swab	Cassette	20T
	H138-21U			Swab/Urine	Cassette	25T
	H139-11S	Chlamydia Trachomatis Ag Rapid Test	CE ₀₁₂₃ NMPA	Swab	Strip	50T
	H139-21S		CE ₀₁₂₃ NMPA	Swab	Cassette	20T
	H139-21U	Chlamydia Trachomatis Ag Test		Swab/Urine	Cassette	25T
	H1310-41S	Gonorrhoeae/Chlamydia Trachomatis Ag Rapid Test		Swab	Panel	20T
	H1311-11W	Chlamydia Trachomatis IgG Ab Rapid Test		WB/S/P	Strip	50T
	H1311-21W			WB/S/P	Cassette	20T
Respiratory Tests	H191-21S	Influenza A/B Ag Rapid Test (Tri-line)	CE NMPA	Swab	Cassette	25T
	H192-21S	Strep A Ag Rapid Test	CE	Swab	Cassette	25T
	H193-21S	ADV Respiratory Adenovirus Ag Rapid Test	CE	Swab	Cassette	25T
	H194-21S	RSV Respiratory Syncytial Virus Ag Rapid Test	CE	Swab	Cassette	25T
	H195-41S	ADV/RSV Ag 2-in-1 Rapid Test	CE	Swab	Panel	20T
	H196-41S	Influenza A/B/RSV/ADV Ag 4-in-1 Rapid Test	CE	Swab	Panel	20T
	H197-21S	MP Mycoplasma Pneumonia Ag Rapid Test	CE	Swab	Cassette	25T
	H198-21WM	MP Mycoplasma Pneumonia IgM Ab Rapid Test	CE NMPA	WB/S/P	Cassette	25T
	H198-21WG	MP Mycoplasma Pneumonia IgG Ab Rapid Test	CE NMPA	WB/S/P	Cassette	25T

	Cat. No.	Product Description	Qualification	Specimen	Format	Kit Size
Respiratory Tests	H198-41W	MP Mycoplasma Pneumonia IgM/IgG Ab Rapid Test	CE	WB/S/P	Panel	25T
	H199-21WM	Cpn Chlamydia Pneumonia IgM Ab Rapid Test	NMPA	WB/S/P	Cassette	25T
	H199-21WG	Cpn Chlamydia Pneumonia IgG Ab Rapid Test	NMPA	WB/S/P	Cassette	25T
	H100A-21W	SARS-CoV-2 IgM/IgG Ab Rapid Test (Tri-line)	CE NMPA	WB/S/P	Cassette	25T
	H100N-21W	SARS-CoV-2 Neutralizing Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H100V-41W	SARS-CoV-2 Ab/Neutralizing Ab Rapid Test (Tri-line)	CE	WB/S/P	Panel	25T
	H100D-21W	SARS-CoV-2 RBD Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H100DG-21W	SARS-CoV-2 RBD IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H100G-21S	SARS-CoV-2 Ag Rapid Test (For professional use)	CE	Swab	Cassette	1T/5T/25T
	H100G	SARS-CoV-2 Ag Rapid Test (For self test)	CE ₀₁₂₃ NMPA	Swab	Cassette	1T/5T/25T
	H100S-21A	SARS-CoV-2 Ag Saliva Rapid Test	CE	Saliva	Cassette	20T
	H100C-41S	SARS-CoV-2/FLU A+B Ag Rapid Test	CE	Swab	Panel	20T
	H107-21S	SARS-CoV-2/Influenza A/B Ag Rapid Test	CE	Swab	Cassette	25T
	H108-41S	SARS-CoV-2/FLU A+B/RSV Ag Rapid Test	CE	Swab	Panel	20T
	H109-41S	SARS-CoV-2/FLU A+B/RSV/ADV Ag Rapid Test	CE	Swab	Panel	20T
	H1010-41S	SARS-CoV-2/FLU A+B/RSV/ADV/MP Ag Rapid Test		Swab	Panel	20T
	H1910-41W	Respiratory 7-in-1 IgM Antibody Rapid Test	NMPA	WB/S/P	Panel	20T
	H1911-21S	hMPV Ag Rapid Test		Swab	Cassette	25T
	H1912-21S	RhV Rhinovirus Ag Rapid Test		Swab	Cassette	25T
	H1914-21S	HPIV Human Parainfluenza Virus Ag Rapid Test		Swab	Cassette	25T
	H1424-11WG	TB Tuberculosis IgG Ab Rapid Test	CE NMPA	WB/S/P	Strip	50T
	H1424-21WG	TB Tuberculosis IgG Ab Rapid Test	CE NMPA	WB/S/P	Cassette	25T
	H1424-21W	TB Tuberculosis IgG/IgM Ab Rapid Test (Tri-line)	CE	WB/S/P	Cassette	25T
	H1437-21W	Measles IgG/IgM Ab Rapid Test		WB/S/P	Cassette	25T
Prenatal and Postnatal Care Tests	H181-21WM	TOXO Toxoplasma IgM Ab Rapid Test	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H181-21WG	TOXO Toxoplasma IgG Ab Rapid Test	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H181-21W	TOXO Toxoplasma Ab IgM/IgG Test (Tri-line)	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H181-41W	TOXO Toxoplasma Ab IgM/IgG Test	CE ₀₁₂₃	WB/S/P	Panel	25T
	H182-21WM	RV Rubella IgM Ab Rapid Test	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H182-21WG	RV Rubella IgG Ab Rapid Test	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H182-21W	RV Rubella Ab IgM/IgG Test (Tri-line)	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H182-41W	RV Rubella Ab IgM/IgG Test	CE ₀₁₂₃	WB/S/P	Panel	25T
	H183-21WM	CMV Cytomegalovirus IgM Ab Rapid Test	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H183-21WG	CMV Cytomegalovirus IgG Ab Rapid Test	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H183-21W	CMV Cytomegalovirus Ab IgM/IgG Test (Tri-line)	CE ₀₁₂₃	WB/S/P	Cassette	25T
	H-183-41W	CMV Cytomegalovirus Ab IgM/IgG Test	CE ₀₁₂₃	WB/S/P	Panel	25T
	H184-21WM	HSV-I Herpes Simplex Virus I IgM Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H184-21WG	HSV-I Herpes Simplex Virus I IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H184-21W	HSV-I Herpes Simplex Virus I Ab IgM/IgG Test (Tri-line)	CE	WB/S/P	Cassette	25T
	H184-41W	HSV-I Herpes Simplex Virus I Ab IgM/IgG Test	CE	WB/S/P	Panel	25T
	H185-21WM	HSV-II Herpes Simplex Virus II IgM Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H185-21WG	HSV-II Herpes Simplex Virus II IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H185-21W	HSV-II Herpes Simplex Virus II Ab IgM/IgG Test (Tri-line)	CE	WB/S/P	Cassette	25T
	H185-41W	HSV-II Herpes Simplex Virus II Ab IgM/IgG Test	CE	WB/S/P	Panel	25T
	H186-41WM	TORCH 5-in-1 IgM Ab Rapid Test	NMPA	WB/S/P	Panel	20T
	H186-41WG	TORCH 5-in-1 IgG Ab Rapid Test	NMPA	WB/S/P	Panel	20T
	H186-41W	TORCH 10-in-1 IgM/IgG Ab Rapid Test		WB/S/P	Panel	20T

Rapid Test

	Cat. No.	Product Description	Qualification	Specimen	Format	Kit Size
Insect-borne Disease Tests	H141-21W	Dengue Ag NS1 Rapid Test	CE	WB/S/P	Cassette	25T
	H142-21W	Dengue IgM/IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H143-41W	Dengue NS1/IgM/IgG Rapid Test	CE	WB/S/P	Panel	20T
	H146-11W	Malaria Pf Ag Rapid Test	CE	Whole Blood	Strip	50T
	H146-21W		CE	Whole Blood	Cassette	25T
	H147-11W	Malaria Pv Ag Rapid Test	CE	Whole Blood	Strip	50T
	H147-21W		CE	Whole Blood	Cassette	25T
	H148-11W	Malaria Pan Ag Rapid Test	CE	Whole Blood	Strip	50T
	H148-21W		CE	Whole Blood	Cassette	25T
	H149-11W	Malaria Pf/Pv Ag Rapid Test	CE	Whole Blood	Strip	50T
	H149-21W		CE	Whole Blood	Cassette	25T
	H1410-11W	Malaria Pf/Pan Ag Rapid Test	CE	Whole Blood	Strip	50T
	H1410-21W		CE	Whole Blood	Cassette	25T
	H1433-41W	Malaria Pf/Pv/Pan Ag Rapid Test		Whole Blood	Cassette	25T
	H1434-41W	Malaria Pf (HRP2/pLDH) Ag Rapid Test		Whole Blood	Cassette	25T
	H1413-21P	Chikungunya IgM Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H1413-21WT	Chikungunya Ab IgM/IgG Rapid Test (Tri-line)		WB/S/P	Cassette	25T
	H1413-21W	Chikungunya Ag Rapid Test		WB/S/P	Cassette	25T
	H1417-21P	Leishmania IgM/IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H1431-21P	Filariasis IgM/IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H1435-21W	Tsutsugamushi IgM Ab Rapid Test		WB/S/P	Cassette	25T
	H1436-21W	Chagas IgG Ab Rapid Test		WB/S/P	Cassette	25T
Enteric Infection Tests	H144-11W	H.Pylori Ab Rapid Test	CE NMPA	WB/S/P	Strip	50T
	H144-21W		CE NMPA	WB/S/P	Cassette	25T
	H144-21WT	H.Pylori IgG/IgM Ab Rapid Test (Tri-line)		WB/S/P	Cassette	25T
	H145-11F	H.Pylori Ag Rapid Test	CE	Feces	Strip	20T
	H145-21F		CE	Feces	Cassette	20T
	H1411-21W	Typhoid IgM/IgG Ab Rapid Test (Tri-line)	CE	WB/S/P	Cassette	25T
	H1412-21F	Typhoid Ag Rapid Test	CE	Feces	Cassette	20T
	H1414-21F	Cholerae O1 Ag Rapid Rapid Test	CE	Feces	Cassette	20T
	H1415-21F	Cholerae O139 Ag Rapid Rapid Test	CE	Feces	Cassette	20T
	H1416-41F	Cholerae O1/O139 Ag Rapid Test		Feces	Panel	20T
	H1419-11F	Adenovirus Ag Rapid Test	CE NMPA	Feces	Strip	20T
	H1419-21F		CE NMPA	Feces	Cassette	20T
	H1420-11F	Rotavirus Ag Rapid Test	CE NMPA	Feces	Strip	20T
	H1420-21F		CE NMPA	Feces	Cassette	20T
	H1421-41F	Adenovirus/Rotavirus Ag Rapid Test	CE	Feces	Panel	20T
	H1425-21F	Norovirus GI/GII Ag Rapid Test (Tri-line)	CE	Feces	Cassette	25T
	H1426-41W	EV71/CA16 IgM Ab Rapid Test	CE	WB/S/P	Panel	25T
	H1427-41F	Rotavirus/Adenovirus/Norovirus GI/Norovirus GII/Astrovirus Rapid Test	CE	Feces	Panel	20T
	H1430-21P	EV71 Enterovirus 71 IgM Ab Rapid Test		WB/S/P	Cassette	25T
	H1432-41F	Rotavirus/Adenovirus/Norovirus GI /Norovirus GII Rapid Test		Feces	Panel	20T
Zoonotic Infection Tests	H1418-21P	Leptospira IgM/IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T
	H1428-21W	Human Brucella IgM/IgG Ab Rapid Test		WB/S/P	Cassette	25T
	H1429-21W	Human Brucella Ag Rapid Test		WB/S/P	Cassette	25T
	H1913-21S	Influenza A Virus (H5N1) Ag Rapid Test		Swab	Cassette	25T
	H200S	Monkeypox Virus Ag Rapid Test	CE	Swab	Cassette	25T
	H200B	Monkeypox Virus IgM/IgG Ab Rapid Test	CE	WB/S/P	Cassette	25T

Drug of Abuse Tests

Cat. No.	Product Description	Qualification	Specimen	Cut-Off	Format	Kit Size
H151-11U	AMP Amphetamine Rapid Test	CE	Urine	1000ng/mL,500ng/mL,300ng/mL	Strip	50T
H151-21U		CE	Urine	1000ng/mL,500ng/mL,300ng/mL	Cassette	25T
H151-21A		CE	Saliva	50ng/mL	Cassette	1T/15T/20T 25T/40T
H152-11U	BAR Barbiturates Rapid Test	CE	Urine	200ng/mL,300ng/mL	Strip	50T
H152-21U		CE	Urine	200ng/mL,300ng/mL	Cassette	25T
H152-21A			Saliva	50ng/mL	Cassette	15T
H153-11U	BUP Buprenorphine Rapid Test	CE	Urine	10ng/mL	Strip	50T
H153-21U		CE	Urine	10ng/mL	Cassette	25T
H153-21A			Saliva	5ng/mL	Cassette	15T
H154-11U	BZO Benzodiazepines Rapid Test	CE	Urine	300ng/mL,200ng/mL	Strip	50T
H154-21U		CE	Urine	300ng/mL,200ng/mL	Cassette	25T
H154-21A		CE	Saliva	20ng/mL	Cassette	1T/15T/20T 25T/40T
H155-11U	COC Cocaine Rapid Test	CE	Urine	300ng/mL,150ng/mL	Strip	50T
H155-21U		CE	Urine	300ng/mL,150ng/mL	Cassette	25T
H155-21A		CE	Saliva	20ng/mL	Cassette	1T/15T/20T 25T/40T
H156-11U	COT Continine Test	CE	Urine	200ng/mL	Strip	50T
H156-21U		CE	Urine	200ng/mL	Cassette	25T
H156-21A			Saliva	20ng/mL	Cassette	15T
H157-11U	EDDP 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine Rapid Test	CE	Urine	100ng/mL,300ng/mL	Strip	50T
H157-21U		CE	Urine	100ng/mL,300ng/mL	Cassette	25T
H157-21A			Saliva	20ng/mL	Cassette	15T
H158-11U	KET Ketamine Rapid Test	CE	Urine	1000ng/mL,500ng/mL	Strip	50T
H158-21U		CE	Urine	1000ng/mL,500ng/mL	Cassette	25T
H158-21A			Saliva	50ng/mL,100ng/mL,300ng/mL	Cassette	15T
H159-11U	MDMA Ecstasy Rapid Test	CE	Urine	1000ng/mL,500ng/mL	Strip	50T
H159-21U		CE	Urine	1000ng/mL,500ng/mL	Cassette	25T
H159-21A		CE	Saliva	50ng/mL	Cassette	1T/15T/20T 25T/40T
H1510-11U	MET Methamphetamine Rapid Test	CE	Urine	1000ng/mL,500ng/mL	Strip	50T
H1510-21U		CE	Urine	1000ng/mL,500ng/mL	Cassette	25T
H1510-21A		CE	Saliva	35ng/mL,50ng/mL	Cassette	1T/15T/20T 25T/40T
H1511-11U	MOP Morphine Rapid Test	CE	Urine	300ng/mL,100ng/mL	Strip	50T
H1511-21U		CE	Urine	300ng/mL,100ng/mL	Cassette	25T
H1511-21A		CE	Saliva	15ng/mL	Cassette	1T/15T/20T 25T/40T
H1512-11U	MTD Methadone Rapid Test	CE	Urine	300ng/mL,200ng/mL	Strip	50T
H1512-21U		CE	Urine	300ng/mL,200ng/mL	Cassette	25T
H1512-21A		CE	Saliva	30ng/mL	Cassette	1T/15T/20T 25T/40T
H1513-11U	OPI Opiates Rapid Test	CE	Urine	2000ng/mL	Strip	50T
H1513-21U		CE	Urine	2000ng/mL	Cassette	25T
H1513-21A		CE	Saliva	40ng/mL	Cassette	1T/15T/20T 25T/40T
H1514-11U	OXY Oxycodone Rapid Test	CE	Urine	300ng/mL,100ng/mL	Strip	50T
H1514-21U		CE	Urine	300ng/mL,100ng/mL	Cassette	25T
H1514-21A			Saliva	20ng/mL	Cassette	15T
H1515-11U	PCP Phencyclidine Rapid Test	CE	Urine	25ng/mL,10ng/mL	Strip	50T
H1515-21U		CE	Urine	25ng/mL,10ng/mL	Cassette	25T
H1515-21A		CE	Saliva	10ng/mL	Cassette	1T/15T/20T 25T/40T
H1516-11U	PPX Propoxyphene Rapid Test	CE	Urine	300ng/mL,100ng/mL	Strip	50T
H1516-21U		CE	Urine	300ng/mL,100ng/mL	Cassette	25T
H1516-21A			Saliva	50ng/mL	Cassette	15T

Rapid Test

Cat. No.	Product Description	Qualification	Specimen	Cut-Off	Format	Kit Size
H1517-11U	K2 Synthetic Marijuana Rapid Test	CE	Urine	50ng/mL	Strip	50T
H1517-21U		CE	Urine	50ng/mL	Cassette	25T
H1517-21A			Saliva	25ng/mL	Cassette	15T
H1518-11U	TCA TricyclicAntidepressants Rapid Test	CE	Urine	1000ng/mL,500ng/mL	Strip	50T
H1518-21U		CE	Urine	1000ng/mL,500ng/mL	Cassette	25T
H1518-21A			Saliva	100ng/mL	Cassette	15T
H1519-11U	THC Marijuana Rapid Test	CE	Urine	20ng/mL,25ng/mL,50ng/mL	Strip	50T
H1519-21U		CE	Urine	20ng/mL,25ng/mL,50ng/mL	Cassette	25T
H1519-21A		CE	Saliva	25ng/mL,50ng/mL	Cassette	1T/15T/20T 25T/40T
H1520-11U	TML Tramadol Rapid Test	CE	Urine	100ng/mL,300ng/mL,200ng/mL	Strip	50T
H1520-21U		CE	Urine	100ng/mL,300ng/mL,200ng/mL	Cassette	25T
H1520-21A			Saliva	50ng/mL	Cassette	15T
H1521-11U	FYL Fentanyl Rapid Test	CE	Urine	10ng/mL,100ng/mL,200ng/mL	Strip	50T
H1521-21U		CE	Urine	10ng/mL,100ng/mL,200ng/mL	Cassette	25T
H1521-21A			Saliva	20ng/mL	Cassette	15T
H1521-71U	FYL Fentanyl Urine Rapid Test	FDA	Urine	1ng/mL	Dipcard	1T/25T
H1522-11U	ETG Ethyl Glucuronide Rapid Test	CE	Urine	1000ng/mL	Strip	50T
H1522-21U		CE	Urine	1000ng/mL	Cassette	25T
H1523-11U	SOMA Carisoprodol Rapid Test		Urine	2000ng/mL	Strip	50T
H1523-21U			Urine	2000ng/mL	Cassette	25T
H1524-11U	PGB Pregabalin Rapid Test		Urine	5000ng/mL	Strip	50T
H1524-21U			Urine	5000ng/mL	Cassette	25T
H1525-11U	6-MAM 6-Monoacetylmorphine Rapid Test		Urine	10ng/mL	Strip	50T
H1525-21U			Urine	10ng/mL	Cassette	25T
H1530-11U	K3 Sythetic Marijuana Rapid Test		Urine	10ng/mL	Strip	50T
H1530-21U			Urine	10ng/mL	Cassette	25T
H1531-11U	MQL Methaqualone Rapid Test		Urine	300ng/mL	Strip	50T
H1531-21U			Urine	300ng/mL	Cassette	25T
H1532-11U	LSD Lysergic Acid Diethylamide Rapid Test		Urine	20ng/mL	Strip	50T
H1532-21U			Urine	20ng/mL	Cassette	25T
H1533-11U	ZOL Zolpidem Rapid Test		Urine	50ng/mL	Strip	50T
H1533-21U			Urine	50ng/mL	Cassette	25T
H1536-11U	GAB Gabapentin Rapid Test		Urine	1000ng/mL	Strip	50T
H1536-21U			Urine	1000ng/mL	Cassette	25T
H1537-11U	MDPV Methylendioxypropyrolone Rapid Test		Urine	500ng/mL, 1000ng/mL	Strip	50T
H1537-21U			Urine	500ng/mL, 1000ng/mL	Cassette	25T
H1538-11U	COD Codeine Rapid Test		Urine	200ng/mL, 300ng/mL	Strip	50T
H1538-21U			Urine	200ng/mL, 300ng/mL	Cassette	25T
H1534-71U	NFYL Norfentanyl Urine Rapid Test	FDA	Urine	5ng/mL	Dipcard	1T/25T
H1535-41U	FYL/NFYL Fentanyl/Norfentanyl Urine Rapid Test	FDA	Urine	FYL1ng/mL, NFYL5ng/mL	Panel	1T/25T
H1526-51U	Multi-Drug Rapid Test (1-14 Para)	CE	Urine	-	Panel	20T
H1526-51A	Multi-Drug Rapid Test		Saliva	-	Panel	20T
H1526-71U	Multi-Drug Test Cassette		Urine	-	Cassette	20T
H1526-71A	Multi-Drug Test Cassette		Saliva	-	Cassette	20T
H1527-61U	Multi-Drug Rapid Test (5-18 Para)	CE FDA	Urine	-	L-Cup	1T
H1528-61U		CE	Urine	-	K-Cup	1T
H1529-61A	Multi-Drug Rapid Test		Saliva	-	Cup	1T

Multi-Drug Test Panel available for free combination of following test items:

AMP, BAR, BUP, BZO, COC, COT, EDDP, ETG, FYL, K2, K3, KET, LSD, 6-MAM, MET, MDMA, MOP/OPI, MQL, MTD, NFYL, OXY, PCP, PGB, PPX, SOMA, TCA, THC, TML, ZOL, GAB, MDPV, COD.

*Adulteration tests available on-demand.

Cardiac Marker Tests

Cat. No.	Product Description	Qualification	Specimen	Cut-Off	Format	Kit Size
H161-11W	cTnI Troponin I Rapid Test	CE NMPA	WB/S/P	1ng/mL	Strip	50T
H161-21W		CE NMPA	WB/S/P	1ng/mL	Cassette	25T
H162-11W	Myo Myoglobin Rapid Test	CE	WB/S/P	80ng/mL	Strip	50T
H162-21W		CE	WB/S/P	80ng/mL	Cassette	25T
H163-11W	CK-MB Creatine kinase isozyme Rapid Test	CE	WB/S/P	5ng/mL	Strip	50T
H163-21W		CE	WB/S/P	5ng/mL	Cassette	25T
H164-41W	cTnI/Myo/CK-MB Rapid Test	CE NMPA	WB/S/P	-	Panel	20T

Tumor Marker Tests

Cat. No.	Product Description	Qualification	Specimen	Cut-Off	Format	Kit Size
H171-11W	AFP Alpha Fetoprotein Rapid Test	CE	WB/S/P	20ng/mL	Strip	50T
H171-21W		CE	WB/S/P	20ng/mL	Cassette	25T
H172-11W	CEA Carcinoma Embryonic Ag Rapid Test	CE	WB/S/P	5ng/mL	Strip	50T
H172-21W		CE	WB/S/P	5ng/mL	Cassette	25T
H173-11F	FOB Fecal Occult Blood Rapid Test	CE NMPA	Feces	50ng/mL	Strip	20T
H173-21F		CE NMPA	Feces	50ng/mL	Cassette	20T
H174-11W	PSA Prostate Specific Ag Rapid Test	CE ₀₁₂₃	WB/S/P	4ng/mL	Strip	50T
H174-21W		CE ₀₁₂₃	WB/S/P	4ng/mL	Cassette	25T
H175-41W	AFP/CEA/PSA Rapid Test		WB/S/P	-	Panel	20T
H176-21F	TF Transferrin Rapid Test	CE	Feces	40ng/mL	Cassette	20T
H177-41F	Hemoglobin (Hb)/Transferrin (TF) Rapid Test	CE NMPA	Feces	-	Panel	20T

Alcohol Tests

Cat. No.	Product Description	Qualification	Specimen	Format	Kit Size
H001-11A	One Step Saliva Alcohol Test	CE	Saliva	Strip	5T/10T/25T/50T
H001-21A		CE	Saliva	Cassette	5T/10T/25T
H001-11U	One Step Urine Alcohol Test	CE	Urine	Strip	5T/10T/25T/50T
H001-21U		CE	Urine	Cassette	5T/10T/25T

Inflammation Tests

Cat. No.	Product Description	Qualification	Specimen	Cut-Off	Format	Kit Size
H002-11P	Serum Cholinesterase Test		Serum	-	Strip	50T
H003-21W	RF Rheumatoid Factor Test	CE	WB/S/P	25IU/mL	Cassette	25T
H004-21W	ASO Anti-Antistreptolysin O Test	CE	WB/S/P	200IU/mL	Cassette	25T
H1422-11W	CRP C-Reactive Protein Semi-Quantitative Rapid Test	CE NMPA	WB/S/P	-	Strip	50T
H1422-21W		CE NMPA	WB/S/P	-	Cassette	25T
H1423-11W	CCP Anti-cyclic Citrullinated Peptide Rapid Test	CE NMPA	WB/S/P	-	Strip	50T
H1423-21W		CE NMPA	WB/S/P	-	Cassette	25T
H006-21F	Calprotectin Rapid Test		Feces		Cassette	25T

Rapid Test @

Other Rapid Tests

Cat. No.	Product Description	Qualification	Specimen	Format	Kit Size
H005-21W	25-Hydroxy Vitamin D Rapid Test	CE NMPA	WB/S/P	Cassette	25T
H007-21S	PEP Pepsin Rapid Test		Saliva	Cassette	25T
H008-21W	ABO and RhD Blood Grouping Test (Solid-phase method)		Whole Blood	Cassette	20T
H009-21U	Micro-albumin Rapid Test		Urine	Cassette	25T



TORCH-IgM 5-in-1 Combo Rapid Test



Drug of Abuse Rapid Tests



HIV 1/2 Ab Rapid Test



HBsAg Rapid Test




HCV Ab Rapid Test



Syphilis (TP) Ab Rapid Test


Targeted Polypeptide Gold Immunochromatography Technology

We independently developed a 200nm targeted polypeptide gold, which features a low coefficient of variation in particle size and strong color development, providing a revolutionary solution for the color development and interpretation of rapid test reagents. By innovatively adopting a directional coupling strategy, we can successfully enhance the efficiency of antigen/antibody labeling.



Biotin-Streptavidin Amplification Technology

Established a fusion technology of targeted polypeptide gold and biotin amplification process to achieve multi-level signal amplification, enabling high-sensitivity detection at a relatively low analyte level.



Rapid Test

DOT Kits

	Cat. No.	Product Description	Qualification	Detection	Specimen	Kit Size
Respiratory Tract Tests	H211	TB (Mycobacterium Tuberculosis) Antibody Test	NMPA	IgG	Serum	20T
	H212	MP (Mycoplasma Pneumoniae) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H213	Cpn (Chlamydia Pneumoniae) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H214	RSV (Respiratory Syncytial Virus) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H215	Adv (Adenovirus) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H216	MV (Measles Virus) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H217	PFV (Parainfluenza Virus) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H218	FLU-A (Influenza Virus A) Antibody Test	NMPA	IgM	Serum	20T
	H219	FLU-B (Influenza Virus B) Antibody Test	NMPA	IgM	Serum	20T
	H2110	CVB (Coxsackievirus B) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H2111	EBV (Epstein-Barr Virus) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
	H2112	HPV B19 (Human Parvovirus B19) Antibody Test	NMPA	IgM/IgG	Serum	20T
Sexually Transmitted Disease Tests	H221	Uu (Urea Plasma Urealyticum) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H222	Ct (Chlamydia Trachomatis) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H223	HPV (Human Papilloma Virus) Antibody Test		IgM/IgG	Serum	20T
Hepatitis Tests	H231	HCV Antibody Test		IgG	Serum	20T
	H232	HAV Antibody Test		IgM	Serum	20T
Infectious Disease Tests	H241	Hantavirus Antibody Test	NMPA	IgM/IgG	Serum	20T
	H242	Enterovirus 71 Antibody Test	NMPA	IgM	Serum	20T
	H243	CA16 (Coxsackievirus A16) Antibody Test	NMPA	IgM	Serum	20T
Rheumatic Diseases Tests	H251	ds-DNA (Double-stranded DNA) Antibody Test	NMPA	IgG	Serum	20T
	H252	ss-DNA (Single-stranded DNA) Antibody Test	NMPA	IgG	Serum	20T
	H253	RF (Rheumatoid Factor) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H254	ASO (Anti-streptolysin O) Antibody Test	NMPA	IgG	Serum	20T
	H255	ANA (Antinuclear) Antibody Test	NMPA	IgG	Serum	20T
	H256	ENA (Extractable Nuclear Antigen) Antibody Test	NMPA	IgG	Serum	20T
	H257	CCP (Cyclic citrullinated peptide) Antibody Test	NMPA	IgG	Serum	20T
Immunological Infertility Tests	H261	ASA (Anti-sperm) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
	H262	EM (Anti-endometrial) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
	H263	ACA (Anticardiolipin) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
	H264	AOA (Anti-ovary) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
	H265	HCG (Anti-Human Chorionic Gonadotropin) Antibody Test	NMPA	IgG	Serum	20T
	H266	TA (Anti-Trophoblast Membrane) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
	H267	ZP (Anti-Zona Pellucida) Antibody Test	NMPA	IgA/IgM/IgG	Serum	20T
Other Tests	H271	HP (Helicobacter Pylori) Antibody Test	NMPA	IgM/IgG	Serum	20T
	H272	Allergen Specific IgE Antibody Test kit (Western blotting)		IgE	Serum/Plasma Whole Blood	1T/5T 10T
	H273	Food Specific IgG Antibody Test Kit (Western blotting)	NMPA	IgG	Serum/Plasma Whole Blood	1T/5T 10T



Exclusive Product

Allergen Specific IgE Antibody Test kit (Western blotting)

Identification of Product	Test Item
Comprehensive Group	Peanuts, soy, milk, crab, shrimp, eggs, beef, codfish, wheat, lamb, house dust mites, dust mites, cat dander, dog dander, house dust, cockroaches, Alternaria, willow, ragweed, mugwort
Food Group	Peanuts, soybeans, milk, crab, shrimp, eggs, beef, cod, wheat, lamb
Inhalation Group	House dust mites, dust mites, cat dander, dog dander, house dust, cockroaches, Alternaria, willow, ragweed, mugwort

Food Specific IgG Antibody Test Kit (Western blotting)

Identification of Product	Test Item
Food Group 1	Beef, lamb, peach, corn, soybeans, chicken, goat's milk, milk, peanuts, eggs, mushrooms, crab, oysters, rice, tomatoes, shrimp, pork, wheat, cod, celery
Food Group 2	Beef, peaches, corn, soybeans, chicken, milk, eggs, mushrooms, crabs, rice, tomatoes, shrimp, pork, wheat, cod
Food Group 3	Beef, soy, chicken, milk, eggs, crab, rice, shrimp, wheat, cod

* Allergen Test Panel available for free combination of given test items

Allergen Specific IgE Ab Test Kit (Western blotting)



- **Micro Sample Volume**
- **Instant Result in 10 min**
- **Multiple Sample Coverage**
- **Detection of 20+ Common Allergens**



Food Specific IgG Ab Test Kit (Western blotting)





Contact HIGHTOP Now

Qingdao Hightop Biotech Co., Ltd. No.369 Hedong Road, Hi-tech Industrial
Development Zone, 266112, Qingdao, Shandong, PEOPLE'S REPUBLIC OF CHINA

☎ 0086-532-58710705 ✉ sales@hightopbio.com 🌐 www.hightopbio.com

