

Urinalysis Analyzer

HT-111A

Urinalysis Tests

- FAST ● Sensitivity
- Accuracy



Product Characteristics

- **High Reliability**
 - High brightness cold light source, long service life
- **Simple Operation**
 - User-friendly interface, easily understand and operate
 - Functional partition is scientific and reasonable, easy to manager user information
- **Multiple application scenarios, suitable for small to medium labs**
- **Light & Sensitivity**
 - Small size, saves space
 - New optical sensor, highly sensitive, strong anti-jamming
- **Easy Maintenance**
 - Automatically eliminate fault according to alarms
 - Boot self-test to ensure accurate detection



Qingdao Hightop Biotech Co., Ltd.

HT-111A Urinalysis Analyzer

Project	Technical index
Test Items	LEU, NIT, URO, PRO, pH, BLD, SG, KET, BIL, GLU, VC, Ca, CRE, MCA, A:C
Applicable test paper	8A, 10A, 11A, 12A, 13A, 14A, 4S, 10S
Principle of measurement	Adopts latest optoelectronic technology and microprocessor technology for clinical urine test
Test speed	60/120 reports/hour
Input mode	Keyboard input
Report storage	Regular report: 5,000 QC report: 1,000
External interface	192*64 dots screen
Working pressure	100~240V,50/60Hz
Size	285mm×237mm×171mm
Weight	1.9kg
Note	The urine analyzer is a closed system, which is used with HIGHTOP urinalysis test strips

Urine analyzer strip

Product advantage

- ◎ Complete items, detection items 4-14 optional
- ◎ Imported raw materials, high color intensity, good differentiation, easy to interpret the result
- ◎ Good stability, can be stored at room temperature for 24 months

HT-L Urine analyzer strip

Test item	LEU、 NIT、 URO、 PRO、 pH、 BLD、 SG、 KET、 BIL、 GLU、 VC、 Ca、 CRE、 MCA、 A:C
Applicable instrument	HT-111A
Product specification	50Tests/Box、 100Tests/Box
Validity period	24 Months

Uncut sheet available of any combination of following parameters

compatible QC reagent available on-demand

Contact HIGHTOP Now

Qingdao Hightop Biotech Co., Ltd. No.369 Hedong Road, Hi-tech Industrial Development Zone, 266112, Qingdao, Shandong, PEOPLE'S REPUBLIC OF CHINA

☎ 0086-532-58710705

✉ sales@hightopbio.com

🌐 www.hightopbio.com